

Development of a Model for Functional Studies of ABCG2 (Breast Cancer Resistance Protein) Efflux Employing a Standard BeWo Clone (B24)

Andrew Crowe¹ and Jeffrey A. Keelan²

¹School of Pharmacy and Curtin Health Innovation Research Institute, Curtin University, Perth, Western Australia, Australia.

²School of Women's and Infants Health, University of Western Australia, Crawley, Western Australia, Australia.

ABSTRACT

Human choriocarcinoma-derived BeWo cells express high levels of breast cancer resistance protein (BCRP/ABCG2) with no functional P-glycoprotein (P-gp) (ABCB1) activity, making them a potential model to study bidirectional ABCG2-mediated drug transport. However, the original BeWo clone (B24) available to researchers does not form confluent monolayers with tight junctions required by the model. Our aim was to adapt culture conditions to attempt to generate confluent BeWo monolayers for drug transport studies using the standard B24 clone. BeWo cells (B24; American Type Culture collection [ATCC]) were cultured in six-well plates or polycarbonate millicell inserts in a number of media formulations, growth supplements, and basement membrane substitutes. Cells were examined for confluence by microscopy, and transepithelial electrical resistance (TEER) was measured daily; monolayer permeability was assessed when TEER had stabilized. Optimal growth rates were achieved in culture conditions consisting of Medium 199 (M199) supplemented with epidermal growth factor (EGF; 20 ng/mL), vitamin supplements, and 10% fetal calf serum (FCS) with collagen coating. A TEER of 170 Ω in 0.6 cm² inserts was achieved 2 weeks after seeding under optimal conditions. The cell-impermeable diffusion marker 5(6) carboxy-2,7-dichlorodihydrofluorescein (C-DCDHF) had a permeability coefficient of 3.5×10^{-6} cm/s, indicative of minimal paracellular permeability. ABCG2 expression, as determined by immunoblotting, remained unaffected by confluency. In conclusion, we describe culture conditions for the B24 BeWo clone that facilitate the formation of monolayers with tighter junctions and reduced paracellular transport compared to previously published models. These growth conditions provide a good model of ABCG2-mediated drug transport in a human placental cell line.

INTRODUCTION

Breast cancer resistance protein (BCRP/ABCG2), a member of the ABC family of active efflux proteins, is an important part of the xenobiotic defense system of the body. Along with P-glycoprotein (P-gp)/ABCB1, these two efflux transporters are usually collocated at the same barrier surfaces of the body, such as the blood–brain barrier, kidney, gastrointestinal tract, breast duct, placental syncytium, and fetal capillary endothelium.^{1,2} BCRP has a significant role in preventing xenobiotics and orally active drugs from being absorbed, and from altering the pharmacokinetics of drug disposition, due to its abundance in the gastrointestinal tract.³ Unlike other efflux transporters though, BCRP is a half transporter that requires dimerization of two copies in the plasma membrane for the generation of a functional efflux protein.⁴ While P-gp has a wide array of drug groups that are transported,⁵ BCRP is currently thought to have a smaller subset of substrates. Nevertheless, some important antineoplastic drugs are known to be transported by both BCRP and P-gp, such as daunorubicin and doxorubicin,⁶ while other drugs not affected by P-gp, such as mitoxantrone, ciprofloxacin, bisantrene, and topotecan, are transported by BCRP.^{4,7}

It has proved difficult to determine the characteristics of P-gp or BCRP drug substrates using molecular modeling.⁸ Therefore, experimental characterization of drug efflux properties requires the development and application of *in vitro* models using cell lines with specific efflux properties.⁹ The Caco-2 cell line, for example, expresses P-gp on its apical surface and forms a polarized epithelial monolayer in a culture with tight junctions to limit paracellular transport, making it ideal for bidirectional transport studies of P-gp efflux.^{8,9} Currently, the only equivalent cell line for analysis of BCRP substrates is a Madin Darby canine kidney (MDCK) cell line transfected with human BCRP,¹⁰ which is not widely available.

The placenta is the tissue with the greatest amount of BCRP expression, and placental trophoblast cells, which form the main epithelial barrier between the maternal and fetal circulations, are characterized by high levels of BCRP expression. The BeWo cell line

ABBREVIATIONS: Ap, apical; ATCC, American Type Culture collection; Bas, basolateral; BCA, bicinechonic acid; BCRP, breast cancer resistance protein; cAMP, cyclic adenosine monophosphate; C-DCDHF, 5(6) carboxy-2,7-dichlorodihydrofluorescein; dbcAMP, dibutyl cyclic AMP; DMEM, Dulbecco's modified Eagle's medium; EGF, epidermal growth factor; FCS, fetal calf serum; HBSS, Hank's balanced salt solution; IgG, immunoglobulin G; M199, Medium 199; MDCK, Madin Darby canine kidney; MEM, minimum essential medium; NEAA, nonessential amino acids; PAGE, polyacrylamide gel electrophoresis; PBS, phosphate-buffered saline; P-gp, P-glycoprotein; rhEGF, recombinant human epidermal growth factor; RPMI, Roswell Park Memorial Institute; SDS, sodium dodecyl sulfate; TEER, transepithelial electrical resistance.

is derived from human trophoblast cells and possesses two key trophoblast characteristics: high endogenous expression of BCRP^{11,12} and propensity for cell–cell fusion and differentiation—a process known as syncytialization by which a multicellular epithelial membrane is formed across which all substances must pass in order to leave or enter the fetus.¹³ In addition to expressing significant amounts of BCRP, BeWo cells express little ABCB1/P-gp, making them ideal for BCRP transport studies.^{14,15} Other transporters are present on BeWo cells such as ABCB2 and some of the ABCC series of transporters, but these are not considered as important for drug efflux as BCRP and P-gp. Like many other human cell lines, BeWo cells will grow in most standard culture media in the presence of fetal calf serum (FCS) and are capable of secreting their own basement membrane to adhere to the surface of plastic culture dishes. However, under normal culture conditions these cells do not form confluent cell layers; instead, they aggregate and form discrete colonies.¹⁴ At high seeding concentrations the cell colonies syncytialize and retract from other colonies, creating large gaps that preclude the development of intact monolayers.^{15,16} This characteristic renders BeWo cells unusable for bidirectional transport studies. Although a subclone of BeWo cells (B30) has been shown to form tight junctions and has been used in transport studies, there has been very limited use of this subclone and the tight junctions described still allowed substantial paracellular transport.^{17,18} Thus, our aim was to assess whether it is possible to adapt the culture conditions to encourage the B24 clone of BeWo cells available through national cell banks to display tight junctions and form a monolayer barrier in culture, with the ultimate goal of developing a useable *in vitro* model for BCRP-mediated transport studies.

MATERIALS AND METHODS

Topotecan hydrochloride, 5(6)-carboxy-2,7-dichlorodihydrofluorescein (C-DCDHF), bicinchoninic acid (BCA), collagen type I from calf skin, anti- β -actin monoclonal antibody, the Dulbecco's modified Eagle's medium (DMEM), and Medium 199 (M199) were supplied by Sigma-Aldrich. Fumitremorgin C and forskolin were supplied by Enzo Life Sciences, Inc. PSC-833 was kindly donated by Novartis Pharmaceuticals. The Roswell Park Memorial Institute (RPMI) medium was from Hyclone. GIBCO minimum essential medium (MEM) nonessential amino acids (NEAA; Cat. No. 11140-050), GIBCO MEM vitamin supplement (100 \times ; Cat. No. 11120-052), recombinant human epidermal growth factor (rhEGF) Ham's F-12 Nutrient Mixture with Kaighn's Modification (F12K medium), FCS (Cat. No. 10100-139 and 10100-147), GlutaMAXTM, and peroxidase-labeled goat anti-mouse immunoglobulin G (IgG) antibody were supplied by Life Technologies Australia Pty Ltd. FCS was also supplied by Australian Commonwealth Serum Laboratories. Mouse monoclonal anti-human cytokeratin 7 antibody was from DAKO Australia. Precast electrophoresis gels were from PAGEgel Inc. Anti-P-gp MDR1 (G-1) mouse monoclonal, as well as anti-MRP1 mouse monoclonal (QCRL-1), and anti-MRP2 (M2 III-6) mouse monoclonal antibodies were supplied by SantaCruz Biotechnology, while anti-BCRP monoclonal antibody (BXP-21) was purchased from Chemicon (Life Technologies Australia Pty Ltd.).

SuperSignal West-Pico Substrate was from Pierce Chemical. All other reagents were supplied by Sigma-Aldrich Pty. Ltd.

Immunohistochemistry

BeWo cells growing on 6- and 24-well plates or 0.6-cm² inserts were visualized by staining with anti-cytokeratin 7 to determine levels of confluence. Cells were fixed in freshly prepared 4% paraformaldehyde in phosphate-buffered saline (PBS) for 10 min and subsequently blocked with 5% casein in PBS for 1 h. After incubation with primary antibody for 1 h, cells were washed in PBS and peroxidase-labeled secondary anti-mouse IgG antibody was added. After an hour, cells were subject to four stringent washes with PBS and color was developed using the 3,3'-diaminobenzide peroxidase substrate kit (Vector Labs); nickel was omitted for the staining of cells on the 6- and 24-well plates, but included when cells cultured on inserts were stained. Images were collected on a Nikon DS-5Mc camera attached to a Nikon Ti inverted phase microscope using NIS imaging software (Nikon Corp.).

Protein Determination

BeWo cells were scraped and homogenized in a protein lysis buffer consisting of 120 mM NaCl, 24 mM Tris base, 1% Triton-X-100, 0.1% sodium dodecyl sulfate (SDS), and 0.2% deoxycholic acid in water, pH 7.4. The lysis buffer was supplemented with Sigma Ultrafast general protein inhibitor tablets just before use. Protein concentrations in the lysates were determined using a micro-BCA method adapted for use with multiwell plates on a VersaMax plate reader (Molecular Devices). Curve fitting and data extrapolation were performed using SoftMax Pro40 software (Molecular Devices).

Western Blotting for Efflux Protein Detection

Thirty micrograms of proteins was added to each well of a 4%–12% Tris-glycine gel and subjected to sodium dodecyl sulfate–polyacrylamide gel electrophoresis (SDS-PAGE) before transfer onto nitrocellulose membranes. They were washed in Tris-buffered saline/0.05% Tween 20, and then probed with monoclonal primary antibodies directed against P-gp MDR1 (G-1), BCRP (BXP-21), MRP1 (QCRL-1), and MRP2 (M2 III-6) diluted to a final concentration of 0.5 μ g/mL. In parallel, blots were also probed with mouse anti- β -actin antibody at 1:10,000 dilution incubated at room temperature for 1.5 h with gentle rocking. After addition of horseradish peroxidase-labeled goat anti-mouse IgG (1:5,000 dilution), blots were developed by enhanced chemiluminescence and visualized after 20 min exposure using an ImageQuant 350 image analyzer (GE Healthcare).

Optimization of Culture Conditions

BeWo cells were seeded at 30,000 cells/cm² on 6- or 24-well plates (BD Falcon, BD Biosciences) and at 100,000–125,000 cells/cm² on 0.6 cm² Millicell polycarbonate inserts (Merck Millipore). A variety of culture parameters were then evaluated: (1) use of different media formulations (DMEM, F12K, RPMI, or M199), (2) addition of NEAA and/or MEM vitamin solution, and (3) different batches and suppliers of FCS (CSL, Life Technologies). These combinations were also tested

with and without culture surface coating (collagen, poly-L-lysine, or bovine gelatine). For coating with collagen, tissue culture wells and 0.6-cm² polycarbonate inserts were incubated with 150 μ L collagen I (1 mg/mL) from calf skin dissolved in 0.9 M acetic acid for a minimum of 2 h. The acetic acid solution was aspirated and wells were left in the culture cabinet under ultraviolet (UV) illumination for 30 min and then left to dry overnight. The wells were subsequently washed three times in sterile PBS to remove traces of acetic acid before seeding with cells. To coat with poly-L-lysine or gelatine, plates and inserts were coated in 200 μ L of a 1% solution of either poly-L-lysine or gelatine for 1 h before aspirating, UV treating for 30 min, and subsequently leaving to dry overnight before conducting three washing steps to remove traces of nonbound poly-L-lysine or gelatine. To evaluate the effects of differentiation and syncytialization on growth characteristics such as aggregation and confluence, dibutyl cyclic AMP (dbcAMP; 100 μ M), forskolin (10 μ M), or human epidermal growth factor (10–20 ng/mL) were added 24 or 48 h after initial seeding on the wells. Cells were incubated for up to 12 days, and were inspected daily for changes in morphology and confluence. Trans-epithelial electrical resistance (TEER) across the inner and outer compartments of the insert chambers was measured daily using an EVOM meter and the ENDOHM 12 chamber (World Precision Instruments).

Bidirectional Transport

Filter inserts were transferred to fresh 24-well plates for the transport studies. The studies were conducted using an assay medium consisting of Hank's balanced salt solution (HBSS) supplemented with both glucose (Ajax chemicals) and HEPES to give final concentrations of 20 and 10 mM, respectively. The pH was adjusted to 7.4 using 1 M NaOH.

Cells were incubated in a prewarmed assay medium with or without an efflux inhibitor for 30 min at the correct pH, and then rinsed in the same medium. TEER was measured immediately before the start of the experiment and the assay medium with or without inhibitors were placed in the receiver chambers. We used topotecan hydrochloride (20 μ M), a typical BCRP substrate, rhodamine 123 (5 μ M), a typical P-gp substrate, or C-DCDHF, a paracellular transported compound, added to the donor chamber of each well (apical chamber [Ap] for apical to basolateral studies, and basolateral chamber [Bas] for basolateral to apical studies). The Ap and Bas chambers received 0.3 and 0.6 mL of medium, respectively. Samples were removed from the receiver chamber at various time points over a 3-h period. Constant volumes were maintained by adding the prewarmed medium to the receiver chambers to maintain an equilibrium pressure differential between the volumes in the donor and receiver chambers. The selective BCRP inhibitor fumitremorgin C (5 μ M) was used to block functional BCRP activity, while the P-gp inhibitor PSC-833 (4 μ M) was also used to provide comparative non-BCRP efflux inhibition data. During the transport study these modifying agents were also present in the donor chamber with our test compounds and in the receiver chamber, at the same concentrations as stated above. Fluorescence was measured in 96-well plates using a

FLUOstar Optima fluorescence plate reader. Filters used included a 485-nm excitation filter and a 520-nm emission filter for C-DCDHF and rhodamine 123, while a 355-nm excitation filter and a 520-nm emission filter were used for topotecan detection similar to that used previously (Table 1).¹⁹

Data Analysis

Drug transport through cell monolayers was calculated both as a simple amount passing the monolayer per min, which would vary depending on the concentration used in the donor compartment, and as an apparent permeability coefficient as calculated previously.²⁰ Briefly, our permeability coefficient equation (derived from Eqs. 1 and 2) allows for changes in the concentration of test compound in the donor chamber as the permeability experiment progressed.

$$\text{Cl vol} = \frac{A_a}{\left\{ \frac{[C_{d_0} \times V_d - [A_a + (C_a \times V_s)] - (A_c) \times n / n_{\text{fin}}]}{V_d} \right\}} \quad (1)$$

$$P_{\text{eff}} = \frac{\text{Cl vol} / n}{A} \quad (2)$$

where Cl. vol is clearance volume (mL); n , time (min); n_{fin} , final time point (min); A_a , amount in acceptor compartment at time n (pmol); C_{d_0} , concentration in donor compartment at time 0 (nM); V_d , volume of donor compartment (mL); C_a , concentration in acceptor compartment at previous sample point (nM), V_s , sample volume of previous time point (mL); A_c , amount of compound associated with cells after cells are removed from filters at the end of the study (pmol); A , area of monolayer (cm²); P_{eff} , effective permeability (cm/s).

Results in this study are presented as the mean \pm standard error of the mean (SEM). Significant differences between values were examined using the Student's two-tailed unpaired t -test with results considered significant if $P < 0.05$.

RESULTS

Optimization Studies

The culture of BeWo cells over 1 week in the various medium formulations and supplements revealed marked differences in growth. The DMEM was the poorest performing medium with very little growth observed; the RPMI medium was marginally better, but neither resulted in sustained cell growth past 3 days, leaving more than 50% of the available surface unoccupied by day 5 (results not shown). Cells grown in the F12K and M199 medium resulted in the formation of isolated clusters of cells, achieving \sim 70% coverage. NEAA supplementation of the F12K medium increased cell growth marginally, while MEM vitamin supplementation also enhanced attachment and growth. No differences were observed between the different batches of FCS used. Subsequent studies employed F12K media with NEAA/vitamin supplementation or the M199 medium with added vitamins.

With respect to coating options, collagen coating promoted denser cell clusters and enhanced confluency (Fig. 1). Gelatin was less effective than collagen, and poly-L-lysine was completely ineffective. Treatment with forskolin or dbcAMP, which promotes

Table 1. Workflow for Optimized BeWo Bidirectional Transport Study

Step	Parameter	Value	Description
1	Collagen coat inserts	150 μ L	1 mg/mL collagen I
2	Incubate	2 h	Room temp
3	Wash	300 μ L	PBS
4	Wait	overnight	In hood under ultraviolet light
5	Prepare growth medium	200 mL	M199 + vitamins + FCS + Glutamax
6	BeWo cell seeding	125,000 cells/mL	Trypsinize cells
7	Wait	24 h	Attachment phase
8	EGF medium	5 μ L	rhEGF in growth medium
9	Feed cells	12–14 days	Replace medium every 48 h
10	TEER	170 Ω	EndOHM probe
11	BCRP inhibitors	5 μ M	Fumitremorgin C
12	BCRP substrate assay	30–180 min	Use of 20 μ M topotecan hydrochloride
13	Substrate detection	355 nm/520 nm (ex/em)	96-well fluorescent plate reader

Step Notes

- 1 mg/mL collagen I dissolved in 0.9 M acetic acid spread over 0.6-cm² millicel polycarbonate inserts in 24-well plates.
- Use Class II biological hood under ultraviolet light, aspirate all solution in inserts.
- Wash off any acetic acid present around the inserts at least three times.
- Wait overnight before using inserts. Allow to air-dry in sterile hood under ultraviolet light.
- Generate final concentration of 10% FCS, 2 mM glutamax, and 2 mL of a 100 \times minimum essential medium vitamin mix solution.
- Trypsinize cells for 5 min. Need cells from 75-cm² flask to generate cell density to cover a 24-well study.
- Allow full day for cells to plate down in humid 37°C incubator with 5% CO₂.
- Add rhEGF at 20 ng/mL from 100 μ g/mL stocks to the medium immediately before feeding cells.
- Maintain cells in 37°C incubator with 5% CO₂.
- Measure TEER using an EVOM meter with an insert probe (EndOHM 12).
- Replace the medium with HBSS and incubate some wells with the BCRP inhibitor.
- Add topotecan in HBSS with or without BCRP inhibitors to donor sides and collect aliquots from receiver side for up to 3 h.
- Apparent permeabilities (P_{app}) can be determined for each direction of transport from apical to basolateral and basolateral to apical.

BCRP, breast cancer resistance protein; EGF, epidermal growth factor; FCS, fetal calf serum; HBSS, Hank's balanced salt solution; M199, Medium 199; PBS, phosphate-buffered saline; rhEGF, recombinant human epidermal growth factor; TEER, transepithelial electrical resistance.

syncytialization, resulted in the cell aggregates condensing and retracting from each, reducing confluency, and increasing the spaces between the cells (Fig. 1). On the other hand, rhEGF (10 ng/mL) was effective at increasing the growth rate without inducing aggregation or retraction (Fig. 1). The combination of the M199 medium supplemented with rhEGF, vitamins, and collagen coating achieved >99% surface coverage. Increasing the rhEGF concentration to 20 ng/mL resulted in complete coverage with no visible spaces evident. Under these conditions, after 1 week of culture, domes of BeWo cells were evident (Fig. 2), which continued to expand every day

postconfluence. This suggested that the BeWo cells had formed sheets of tightly joined cells capable of generating hydrostatic pressure via differences between apical versus basal secretion. Dome formation is thought to occur from fluid being pumped underneath the cells causing them to lift off from the base of culture plates. In permeability studies, where cells are cultured on a suspended membrane insert, there is an outlet for fluid to move through the cell and out into the Bas chamber, preventing dome formation on cells grown on filter inserts. Thus, although dome formation is an interesting phenomenon on cells grown on the base of plastic plates, these events do not occur on membrane filter grown cells.

Permeability Studies

Cells were seeded on inserts at 100,000–125,000 cells/cm², and TEER was measured daily for 12 days to assess the effects of culture optimization on cellular membrane permeability. Cell cultured at 100,000 cells/cm² exhibited variable TEER measured between replicates of 30%. Therefore, a higher cell density of 125,000 cells/cm² was chosen, resulting in more consistent tight junctions between replicates. Cells grown in the M199 medium plus epidermal growth factor (EGF)/vitamins on collagen-coated inserts clearly exhibited much higher TEER values than the other growth conditions (Fig. 3A), with readings above 165 Ω recorded. After 12 days, TEER values plateaued (data not shown), suggesting that optimal tight junction formation had been reached. Transport of the hydrophilic paracellular transport marker C-DCDHF (23 μ M) was assessed on day 12 of culture. Ap-to-Bas transport rates were 18.3 ng/(cm²·min) (66 \times 10⁻⁶ cm/s) in nonseeded collagen-coated inserts (*i.e.*, no BeWo cells), 3.2 ng/(cm²·min) (6.0 \times 10⁻⁶ cm/s) with BeWo cells cultured on collagen-coated inserts in the F12K medium plus EGF, and 1.9 ng/(cm²·min) (3.5 \times 10⁻⁶ cm/s) with cells grown in the M199 medium plus EGF (Fig. 3B, C).

Topotecan hydrochloride is a known BCRP substrate. Bidirectional transport studies through confluent BeWo cells resulted in Ap to Bas transport of 6.0 \times 10⁻⁶ cm/s (2.77 ng/[cm²·min]), while Bas to Ap transport was 7.8 \times 10⁻⁶ cm/s (3.67 ng/[cm²·min]) (Fig. 4). Each transport analysis was done in triplicate (Fig. 4); error bars representing SEM values are shown, but due to their size are obscured






Fig. 1. Microscopic analysis of BeWo cell growth after 7 days in varying culture conditions. **(A)** The F12K medium, **(B)** M199. BeWo cultures were fixed and stained with anti-human cytokeratin antibody with immunoperoxidase detection as described in Materials and Methods. Cells were grown on tissue culture grade plastic (**a–d**) or collagen 1-coated plastic (**e–h**) in control conditions (**a, e**) or in the presence of 10 μ M forskolin (**b, f**), 100 μ M dbcAMP (**c, g**), or 10 ng/mL EGF (**e, h**) from day 2 until day 7 of culture. Images, shown at 100 \times magnification, were taken on a Nikon Ti inverted microscope. EGF, epidermal growth factor; dbcAMP, dibutyl cyclic AMP; M199, Medium 199.



and difficult to see. The difference in the rate of Ap versus Bas to-potecan transport was significant with $P < 0.005$. Addition of fumitremorgen C, a potent BCRP inhibitor, abolished the difference in transport direction ($P < 0.773$) (Fig. 4). The P-gp inhibitor PSC-833 was also used, but Bas to Ap transport remained significantly higher than Ap to Bas transport ($P < 0.024$). To illustrate specificity for BCRP transport, additional studies using the P-gp substrate rhodamine 123 were performed. Transport of rhodamine 123 through BeWo cells layers was shown to be 6.9×10^{-6} cm/s in both direction without the addition of any inhibitors (Fig. 5) ($P < 0.967$). The addition of 4 μ M PSC-833 did not significantly alter transport rates ($P < 0.499$), nor did the use of fumitremorgen C ($P < 0.889$), suggesting that P-gp substrates are not actively transported in this model.

Fig. 2. Dome formation in confluent cultures of BeWo cells grown for 7 days in 6-well plates coated with collagen 1 in M199 supplemented with 10 ng/mL EGF and vitamin supplement. 100 \times magnification on Nikon Ti inverted phase microscope.



Western blotting was performed to identify levels of ABCB1 (P-gp), ABCG2 (BCRP), ABCC1 (MRP1), and ABCC2 (MRP2) proteins. BCRP was readily detected in the BeWo cells used in our study (Fig. 6A). No visible bands were detected for P-gp in the BeWo lysates under any experimental conditions, consistent with previous studies²¹ (Fig. 6B). MRP1 was also not detected, but MRP2 was shown to be expressed, although the strength of the signal was not as strong as that for BCRP. There were no changes in BCRP abundance in the presence or absence of culture factors used to generate confluence and monolayer formation in this cell line.

DISCUSSION

The BeWo cell line has been described in many studies of trophoblast biology, differentiation, and invasion, but seldom used in drug transport studies despite the fact that it expresses a number of functional drug transport proteins, including BCRP.¹⁶ This is almost certainly due to its propensity for aggregation and fusion, which creates large gaps on the culture plate, precluding formation of an intact monolayer capable of effective membrane barrier function.^{17,22} Most investigators employ the B24 BeWo clone for their studies, which is available commercially from public

Fig. 3. (A) Increases in TEER during growth of BeWo cells on collagen-coated 0.6-cm² millicell polycarbonate inserts grown in either the F12K (□, ◇), or the M199 medium (△, ×) supplemented with either 20 ng/mL EGF (◇, △) or 10 μM dbcAMP (□, ×). BeWo cells grown on uncoated inserts in M199 media are also shown (○). The Ap to Bas transport rates (ng/cm²/min) are shown on the right of the graphs (based on the accumulation of C-DCDHF in the basal compartment). **(B)** Comparison of passive permeability of C-DCDHF through confluent BeWo cultures on collagen-coated 0.6-cm² millicell polycarbonate inserts compared to TEER measurements for the inserts. **(C)** Bidirectional transport study of 23 μM C-DCDHF by confluent BeWo cells (day 12) grown in either the F12K medium (□, ◇), or the M199 medium with supplemented EGF (■, ◆). Results are shown in both Ap to Bas (□, ■) and Bas to Ap (◇, ◆) directions for each culture condition. TEER, transepithelial electrical resistance; C-DCDHF, 5(6)-carboxy-2,7-dichlorodihydrofluorescein, Ap, apical; Bas, basolateral.




Fig. 4. Bidirectional transport of topotecan hydrochloride (20 μM) by BeWo cells grown on 0.6-cm² polycarbonate inserts at day 14 of culture. Cells were maintained in simple balanced salt solution for the 3-h study. Fumitremorgin C (5 μM) or vehicle was added 30 min before the start of the experiment and was maintained in both Ap and Bas chambers. Results are shown for topotecan hydrochloride alone (■, ◆) or with fumitremorgin C (□, ◇). Results are shown for both Ap to Bas (□, ■) and Bas to Ap (◇, ◆) directions for each culture condition. Data shown are mean \pm SEM ($n=3$). Error bars may not be visible due to their small size. SEM, standard error of the mean.




Fig. 5. Assessment of bidirectional transport of rhodamine 123 (5 μM) in BeWo cells grown on 0.6-cm² polycarbonate inserts at day 15 of culture. Cells were maintained in simple balanced salt solution for the 3-h study. P-gp inhibitor PSC-833 (4 μM) or vehicle was added 30 min before the start of the experiment and was maintained in both Ap and Bas chambers. Results are shown for rhodamine 123 alone (■, ◆) or with PSC-833 (□, ◇). Results are shown for both Ap to Bas (□, ■) and Bas to Ap (◇, ◆) directions for each culture condition. Data shown are mean \pm SEM ($n=3$). Error bars may not be visible due to their small size. P-gp, P-glycoprotein.




Fig. 6. (A) Western blot analysis of BCRP/ABCG2 protein expression in BeWo cells under different culture conditions. About 30 μg of protein was loaded per lane. Lanes 1–3: Total protein was extracted on day 7 of culture from BeWo cells grown in M199 with and without 20 ng/mL rhEGF for either 24 or 72 h. Lanes 4 and 5: BeWo cells on day 7 of culture grown on collagen-coated wells, with and without rhEGF for 6 days. (B) Western blot analysis of three other active efflux proteins (P-gp, MRP1, and MRP2) in BeWo cells compared to Caco-2, A549, and MCF7 cells used in our laboratory. Lane 1: Prestained MW stds. Lane 2: Caco-2 cells cultured for 21 days. Lane 3: A549 cells cultured for 6 days. Lanes 4 and 5: BeWo cells grown on collagen-coated wells, with and without rhEGF for 6 days. Lane 6: MCF7 cells cultured for 7 days. BCRP, breast cancer resistance protein; rhEGF, recombinant human epidermal growth factor.

cell depositories such as American Type Culture collection (ATCC). Although an early report described (B24) BeWo cell cultures able to produce TEER values of 600 $\Omega \cdot \text{cm}^2$,²³ suggestive of very effective tight junctions, this has not been subsequently replicated and it is widely believed that BeWo cells do not form confluent monolayers. However, one group of researchers, the Schwartz group in Washington, have described a specific BeWo subclone, B30, capable of forming confluency.²² They have conducted several studies using this clone, including studies on iron¹⁷ and linoleic acid transport.¹⁸ Unfortunately, although a few other groups in Europe and the United States have subsequently used this subclone, there is some variability in the effectiveness of tight junction formation from these other groups.^{24,25} Recently, one group from Europe attempted to generate monolayers from the B30 subclone of BeWo cells and were not able to get the tight junctions to increase beyond 14 $\Omega \cdot \text{cm}^2$, with high fluorescein transport rates equivalent to 10% per hour,²⁴ while another study showed tight junctions to be around the 45 $\Omega \cdot \text{cm}^2$ level, and 4% fluorescein transport per hour,²⁵ suggesting that there can be high variability in tight junction formation even with use of a

subclone that has a higher propensity to generate monolayers. In addition, approval to use this subclone must come from specific research labs in the United States, as it is not held in commercial cell banks. In contrast, the B24 clone is widely distributed from official American and European cell line distribution sites, which allows ready worldwide access to laboratories interested in developing a BCRP transport model. Here we have outlined a specific set of culture conditions using standard B24 cells that facilitate the formation of confluent BeWo monolayers with tight junctions that will be useful for the study of trophoblast permeability, substrate specificity, and transport kinetics of efflux proteins such as BCRP/ABCG2 *in vitro*.

We evaluated a variety of culture conditions to arrive at a combination of medium, supplements, and culture surface coating needed to achieve confluency. We optimized the culture medium and vitamin content to achieve optimal growth. We also evaluated the effects of increases in intracellular cyclic adenosine monophosphate (cAMP) levels to promote differentiation and ERK1/2 and p38MAPK phosphorylation^{26–28} based on the reported effects of cAMP on tight junction formation between adjoining BeWo cells.²⁷ However, this strategy merely promoted aggregation and fusion with no benefits on confluence. We also evaluated addition of EGF since this is known to increase cell growth, ERK1/2 and p38MAPK activity, and α 2-integrin expression.²⁹ The addition of EGF proved extremely successful in terms of promoting growth and confluence, particularly at the higher dose (20 ng/mL). It is important to note that EGF was not added until 24 h after seeding to allow normal attachment before growth stimulation, with medium and supplements added every 48 h. In our hands, 12 days of constant exposure to EGF were required to reach peak tight junction formation. Coating with collagen I was also essential for formation of BeWo monolayers in our study. This is despite the findings of Nakatsuji and coworkers that EGF treatment had no impact on the adhesion of BeWo cells to collagen IV coated plates.²⁹ Our observations are unlikely to be related directly to adhesion; it is more likely that cell-to-cell apposition is affected. The differences in EGF response could reflect a difference in the BeWo phenotype when grown on type I and type IV collagen.²⁹

The use of TEER measurements as a marker of barrier formation is a common technique, but subject to considerable laboratory-to-laboratory variance due to lack of consistency and reproducibility in measurement techniques and calibration. The original 1997 study detailing the transport properties of the B30 sub-clone showed fluorescein transport at 10×10^{-6} cm/s and TEER values around $60 \Omega \cdot \text{cm}^2$.¹⁸ In a more recent publication using the B30 clone, Heaton and coworkers reported TEER measurements over $300 \Omega \cdot \text{cm}^2$, but with unexpectedly high rates of transmission of mannitol of 3.5% per hour.¹⁷ Under our optimal growth conditions, BeWo cultures achieved TEER readings of 170Ω , $8.5 \times$ greater than background, while the carboxyfluorescein transport rate was only 1.6% per hour (3.5×10^{-6} cm/s). These TEER inconsistencies may be due to methodological differences, *e.g.*, the use of chopstick electrodes or endohm chambers to measure TEER. We have found that chopstick

electrodes can provide inconsistent and variable data ($> 100 \Omega$ between measures), so now we use fixed chambers that provide more accurate and reproducible results.

Stability of the transport characteristics of the model is an important consideration. In a previous study using the B30 sub-clone on collagen-coated polycarbonate inserts, confluency could not be maintained for more than 60 min.³⁰ In contrast, our model showed linearity in efflux for at least 3 h in BeWo cells maintained in a simple buffered balanced salt solution. Drug transport in polarized epithelia would be expected to exhibit marked differences in directionality. In our model although the A–B versus B–A efflux differential for the BCRP substrate topotecan was significant, it did not exceed 1.3-fold. Criteria established for P-gp-mediated efflux in Caco-2 cells suggest that transport differential should be at least two-fold to achieve a functional model with useful characteristics.³¹ However, this assumes that the extent of polarized expression of the transporters of interest is similar to that of P-gp. The polarity of BCRP expression in BeWo cells has not been assessed, and is currently assumed to be similar to that of the placenta with the majority facing the apical side, although recent findings suggest that BCRP is also present on the maternal side of the placenta.² As our data showed some increased basolateral efflux of topotecan, and use of BCRP inhibitors only decreased Bas to Ap transport. If BCRP was present on both sides of the cells, use of a BCRP inhibitor would have decreased transport in both directions, and this was not observed. In addition, we have also shown that this model does not exhibit P-gp-mediated transport, as no P-gp protein was detected and bidirectional transport studies with a P-gp substrate found no evidence of P-gp-mediated efflux.

Therefore, in this study we have shown that in appropriate culture conditions, standard BeWo cells can generate monolayers with tight junctions for studies of bidirectional transport of high-affinity BCRP substrates and inhibitors. For low-affinity molecules, increases in the polarity of expression or the cellular barrier integrity would be required to provide an adequate level of sensitivity. Measurement of intracellular substrate concentrations might also enhance the utility of the model. While the potential presence of MRP2 must be taken into consideration when using this cell line as a model to study BCRP transport, the absence of P-gp or MRP1 provides a distinct advantage. In addition, the study of other transport proteins and processes, such as the transport of amino acids, fatty acids, and nanoparticles, could also be facilitated by the use of this model.

ACKNOWLEDGMENT

We thank Erin Bolitho from the Pharmacogenetics division within the School of Pharmacy at Curtin University for routine screening of the BeWo cells for mycoplasma and for the preparation of stock concentrations of poly-L-lysine and gelatine solutions.

DISCLOSURE STATEMENT

No competing financial interests exist.

REFERENCES

- Sharom FJ: ABC multidrug transporters: structure, function and role in chemoresistance. *Pharmacogenomics* 2008;9:105–127.
- Keelan JA, Aye IL, Mark PJ, Waddell BJ: ABCA1 and placental cholesterol efflux. *Placenta* 2011;32:708–709.
- Gutmann H, Hruz P, Zimmermann C, Beglinger C, Drewe J: Distribution of breast cancer resistance protein (BCRP/ABCG2) mRNA expression along the human GI tract. *Biochem Pharmacol* 2005;70:695–699.
- Litman T, Brangi M, Hudson E, et al.: The multidrug-resistant phenotype associated with overexpression of the new ABC half-transporter, MXR (ABCG2). *J Cell Sci* 2000;113:2011–2021.
- Ambudkar SV, Kim IW, Sauna ZE: The power of the pump: mechanisms of action of P-glycoprotein (ABCB1). *Eur J Pharm Sci* 2006;27:392–400.
- Ito K, Suzuki H, Horie T, Sugiyama Y: Apical/basolateral surface expression of drug transporters and its role in vectorial drug transport. *Pharm Res* 2005;22:1559–1577.
- Haslam IS, Wright JA, O'Reilly DA, Sherlock DJ, Coleman T, Simmons NL: Intestinal ciprofloxacin efflux; the role of breast cancer resistance protein, BCRP (ABCG2). *Drug Metab Dis* 2011;39:2321–2328.
- Raub TJ: P-glycoprotein recognition of substrates and circumvention through rational drug design. *Mol Pharm* 2006;3:3–25.
- Artursson P, Palm K, Luthman K: Caco-2 monolayers in experimental and theoretical predictions of drug transport. *Adv Drug Deliv Rev* 1996;22:67–84.
- Muenster U, Grieshop B, Ickenroth K, Gnoth MJ: Characterization of substrates and inhibitors for the *in vitro* assessment of Bcrp mediated drug-drug interactions. *Pharm Res* 2008;25:2320–2326.
- Evseenko DA, Murthi P, Paxton JW, et al.: The ABC transporter BCRP/ABCG2 is a placental survival factor, and its expression is reduced in idiopathic human fetal growth restriction. *FASEB J* 2007;21:3592–3605.
- Wang H, Zhou L, Gupta A, et al.: Regulation of BCRP/ABCG2 expression by progesterone and 17 β -estradiol in human placental BeWo cells. *Am J Physiol Endocrinol Metab* 2006;290:E798–E807.
- Drewlo S, Baczynk D, Dunk C, Kingdom J: Fusion assays and models for the trophoblast. *Methods Mol Biol* 2008;475:363–382.
- Masumoto A, Takamoto N, Masuyama H, Akahori Y, Inoue S, Hiramatsu Y: Effects of intermittent high glucose on BeWo choriocarcinoma cells in culture. *J Obstet Gynaecol Res* 2011;37:1365–1375.
- Evseenko DA, Paxton JW, Keelan JA: The xenobiotic transporter ABCG2 plays a novel role in differentiation of trophoblast-like BeWo cells. *Placenta* 2007;28 Suppl A:S116–S120.
- Evseenko DA, Paxton JW, Keelan JA: ABC drug transporter expression and functional activity in trophoblast-like cell lines and differentiating primary trophoblast. *Am J Physiol Regul Integr Comp Physiol* 2006;290:R1357–R1365.
- Heaton SJ, Eady JJ, Parker ML, et al.: The use of BeWo cells as an *in vitro* model for placental iron transport. *Am J Physiol Cell Physiol* 2008;295:C1445–C1453.
- Liu F, Soares MJ, Audus KL: Permeability properties of monolayers of the human trophoblast cell line BeWo. *Am J Physiol* 1997;273(5 Pt 1):C1596–C1604.
- Hubbard KE, Schaiquevich P, Bai F, et al.: Application of a highly specific and sensitive fluorescent HPLC method for topotecan lactone in whole blood. *Biomed Chromatogr BMC* 2009;23:707–713.
- Crowe A, Lemaire M: In vitro and in situ absorption of SDZ-RAD using a human intestinal cell line (Caco-2) and a single pass perfusion model in rats: comparison with rapamycin. *Pharm Res* 1998;15:1666–1672.
- Magnarin M, Rosati A, De Ludicibus S, Bartoli F, Decorti G: Role of ABC transporters in the BeWo trophoblast cell line. *Toxicol Mech Methods* 2008;18:763–769.
- Bode CJ, Jin H, Rytting E, Silverstein PS, Young AM, Audus KL: *In vitro* models for studying trophoblast transcellular transport. *Methods Mol Med* 2006;122:225–239.
- Cerneus DP, van der Ende A: Apical and basolateral transferrin receptors in polarized BeWo cells recycle through separate endosomes. *J Cell Biol* 1991;114:1149–1158.
- Cartwright L, Poulsen MS, Nielsen HM, et al.: *In vitro* placental model optimization for nanoparticle transport studies. *Int J Nanomed* 2012;7:497–510.
- Correia Carreira S, Cartwright L, Mathiesen L, Knudsen LE, Saunders M: Studying placental transfer of highly purified non-dioxin-like PCBs in two models of the placental barrier. *Placenta* 2011;32:283–291.
- Delidakis M, Gu M, Hein A, Vatish M, Grammatopoulos DK: Interplay of cAMP and MAPK pathways in hCG secretion and fusogenic gene expression in a trophoblast cell line. *Mol Cell Endocrinol* 2011;332:213–220.
- Egawa M, Kamata H, Kushiyama A, et al.: Long-term forskolin stimulation induces AMPK activation and thereby enhances tight junction formation in human placental trophoblast BeWo cells. *Placenta* 2008;29:1003–1008.
- Yoshie M, Kaneyama K, Kusama K, et al.: Possible role of the exchange protein directly activated by cyclic AMP (Epac) in the cyclic AMP-dependent functional differentiation and syncytialization of human placental BeWo cells. *Human Reprod* 2010;25:2229–2238.
- Nakatsuji Y, Nishio Y, Tani N, et al.: Epidermal growth factor enhances invasive activity of BeWo choriocarcinoma cells by inducing alpha2 integrin expression. *Endocr J* 2003;50:703–714.
- Mitra P, Audus KL: MRP isoforms and BCRP mediate sulfate conjugate efflux out of BeWo cells. *Int J Pharm* 2010;384:15–23.
- Polli J, Wring S, Humphreys J, et al.: Rational use of *in vitro* P-glycoprotein assays in drug discovery. *J Pharmacol Exp Ther* 2001;299:620–628.

Address correspondence to:

Andrew Crowe, PhD
 BLDG 306, School of Pharmacy
 Curtin University
 Kent St., Bentley
 Perth, WA, Australia

E-mail: a.p.crowe@curtin.edu.au